



Sheet ID PS_07502
Release date : 12/19/2017
Print date : 12/19/2017

Strain name C57BL/6J-Tg(CXCL14)52


Reaction Components

Enzyme TOYOBO KOD FX (Code.No : KFX-101)

Reaction Conditions

[illegible]

Standard PCR



Marker

Gene Ladder 100 (Wako, Cat. 316-06951)

Gel

2% agarose in 1 X TAE (100 V, 30 min)

Comment

PCR products

primer set

product size

product

primer 1	primer 2	650 bp	Tg
primer 3	primer 4	210 bp	Internal control (Tcrd)

Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.