

PCR protocol

RBRC 04050

Strain name C57BL/6N-Tg(CAG-EGFP)14Thg

Primer	primer_name	primer_seq	length	
1	GFP-U279	CGTCCAGGAGCGCACCATCTT	21	Depositor's primer
2	GFP-L624	CTTCTCGTTGGGGTCTTTGCTC	22	Depositor's primer
3	mLC3ex3GT	TGAGCGAGCTCATCAAGATAATCAGGT	27	Mol. Biol. Cell., 15, 1101-1111 (2004).
4	mLC3ex4AG	GTTAGCATTGAGCTGCAAGCGCCGTCT	27	Mol. Biol. Cell., 15, 1101-1111 (2004).


Reaction Components

Enzyme	TaKaRa Ex Taq (Cat. No. RR001A)
--------	---------------------------------

Reaction Conditions

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		11.3			
PCR buffer	10 X	2.0			
dNTP	2.5 mM	1.6			
primer 1	10 pmol/ μ l	1.0			
primer 2	10 pmol/ μ l	1.0			
primer 3	10 pmol/ μ l	1.0			
primer 4	10 pmol/ μ l	1.0			
DNA		1.0			
Taq	5U/ μ l	0.1			
total		20			

step	temp (°C)	time (sec)	note
1	94	180	
2	94	30*	
3	63	30*	
4	72	60*	*Repeat for 30 cycles
5	72	300	



Marker

Gene Ladder 100 (Wako, Cat. 316-06951)

Gel

2% agarose in 1 X TAE (100 V, 30 min)

PCR products

primer_set		product_size	product
primer 1	primer 2	367 bp	Tg (GFP)
primer 3	primer 4	550 bp	Internal control (LC3)

Comment